

Application Note: Acquisition modes

The ADQ series of digitizers supports several modes of acquisition to serve a large variety of measurement situations.

This document is an overview of the acquisition modes. Deviations may occur on specific digitizer models.

1 Introduction

The ADQ series of digitizers supports several acquisition modes. The application note intend to give understanding on how these modes work.

See also the trigger application note for more information about the trigger options. The Data transfer application note gives more information on how to move data from the digitizer to the PC. Also information on how to partition the entire signal chain is in the Data transfer application note.

2 Definitions

#	NAME	DESCRIPTION
1	Record	A sequence of samples
2	Sample	One single data value from the A/D converter
3	ADQ DRAM	DRAM on the digitizer (also SDR14)
4	PC DRAM	RAM in the host PC
5	FPGA RAM	Ram in the FPGA on the digitizer
6	Waveform	See Record

Table 1: Definitions

3 Summary of acquisition modes

REF	NAME	DESCRIPTION	PARAMETERS
4.1	Single-record	One record is collected at the trigger. This is a special case of Multi-record mode	Number of records Record length Pre-trigger Trigger hold-off
4.2	Multi-record	One record is collected at each trigger. The records are stored in ADQ DRAM. Read out is performed after acquisition completed.	Number of records Record length Pre-trigger Trigger hold-off
4.3	Read out while recording	The same as Multi-record, but the read out can start before acquisition is completed.	Number of records Record length Pre-trigger Trigger hold-off
4.4	Continuous multi-record.	The same as Multi-record, but the read out can start at the same time as acquisition. Instead of filling the ADQ DRAM and then read out the data (as for multi-record), the ADQ DRAM is used a circular buffer in a continuous operation.	Record length Pre-trigger Trigger hold-off
4.5	Triggered streaming	At each trigger, a record is collected and sent directly to the host PC.	Record length Pre-trigger Trigger hold-off
4.6	Continuous streaming	At a trigger, a continuous stream of data starts.	Custom

Table 2: Acquisition modes

4 Description of the acquisition modes

4.1 Single-record

This is a special case of Multi-record, [Section 4.2](#), where the number of records is set to 1.

4.2 Multi-record

Multi-record mode is the method for acquiring sets very large records.

Set up a record size and a number of records. The trigger position is controlled via pre-trigger and trigger hold-off parameters. A header is added to each record with trigger information and time stamp. Data is stored in ADQ DRAM and read out is performed when the entire acquisition has completed.

Multi-record mode is preferred when

- Large records, up to millions of samples
- Dense triggers, that is high data rate
- Trigger bursts, that is high average data rate

Limitations to multi-record are

- Total data set must be kept below ADQ DRAM size.




Figure 1: Multi-record.

4.3 Read out while recording¹

Read out while recording is an extension of multi-record that speeds up the read out process.

After setting up a multi-record acquisition, [Section 4.2](#), the read-out is started as soon as the first record has been stored in the ADQ DRAM. This feature puts high requirements on the ADQ DRAM since the ADQ DRAM has to alternate between read and write operations. The real time recording of data from the ADCs is always scheduled with higher priority than the read operation. For high data rate operations, there is not sufficient time for reading data.

The advantage of the read out while recording mode is visible if the triggers are sparse. Then most of the memory content is read out at the end of the batch and the dead time to the next batch can be reduced compared to standard multi-record mode.

Read out while recording multi-record mode is preferred when

- Large records, up to millions of samples
- Sparse triggers, that is the total time for acquiring all records is long.

Limitations to read out while recording multi-record are

- Data rate writing to ADQ DRAM has to leave time for readout.




Figure 2: Multi-record with simultaneous read out.

1. This is sometimes referred to as simultaneous read out and acquire (SAR)

4.4 Continuous Multi-record.

Continuous multi-record enables never-ending acquisition of very large of records.

Continuous multi-record is a further extension of multi-record, [Section 4.2](#), with read-out while acquire, [Section 4.3](#). In the continuous multi-record mode, the trigger rate and the record size is set up to generate a data rate that is lower than the read out speed. The ADQ DRAM is set up to start again from address 0 when full as to provide a cyclic data acquisition. This is a variant of triggered streaming, [Section 4.2](#), which allows very long records.

Continuous multi-record mode is preferred when

- Large records, up to millions of samples
- Long pre-trigger
- Sparse triggers
- The total data set is larger than the ADQ DRAM

Limitations continuous multi-record are

- Average data rate to the ADQ DRAM has to be smaller than the speed of the read out process.




Figure 3: Continuous Multi-record.

4.5 Triggered streaming

Triggered streaming is the method for acquiring very many small records, [Figure 4](#).

Triggered streaming also enables individual triggering of the channels¹, [Figure 5](#).

Triggered streaming do not store data on the ADQ DRAM. All records are recorded in the FPGA FIFO and directly transferred to the host PC. The acquisition is set up with trigger position and record length. A header with time stamp and trigger information is added to each record.

Triggered streaming is preferred when

- Small records, up to kilo samples
- Short re-arm time
- Individual channel trigger and data transfer
- Sparse triggers
- The total data set is larger than the ADQ DRAM

1. Available on ADQ412 only.

Limitations continuous multi-record are

- Average data rate has to be kept below the speed of the read out process.
- Record size has to be smaller than FPGA FIFO




Figure 4: Triggered streaming




Figure 5: Triggered streaming with individual level trigger on ADQ412

4.6 Continuous streaming

Continuous streaming is the method for recording a continuous signal.

At the first trigger (e.g. software command), a continuous stream of data to the host PC is started. The FIFO buffer on the FPGA is limited. One have to make sure that the data rate is reduced to a level at which the host PC can read it. Custom data reduction can be implemented using the ADQ Development Kit. There is no header.

Continuous streaming is preferred when

- There is no records, only streaming data

Limitations continuous multi-record are

- Average data rate has to be limited in the FPGA to match the read out process.




Figure 6: Continuous streaming

5 Summary of available modes (as per 2014-04-14)

ADQ	MULTI RECORD	READ OUT WHILE RECORDING	CONTINUOUS MULTI RECORD	TRIGGERED STREAMING	CONTINUOUS STREAMING
ADQ112	YES	NO	NO	NO	YES
ADQ114	YES	NO	NO	Option	YES
ADQ212	YES	NO	NO	NO	YES
ADQ214	YES	NO	NO	Option	YES
ADQ108	YES	YES	YES	Option	YES
ADQ208	YES	YES	YES	Option	YES
ADQ208D	YES	YES	YES	Option	YES
ADQ1600	YES	YES	YES	Option	YES ¹
ADQ412-1G	YES	YES	YES	YES	YES
ADQ412-3G	YES	YES	YES	YES	YES
ADQ412-4G	YES	YES	YES	YES	YES
SDR14	YES	YES	YES	Option	YES ²

Table 3: Acquisition modes on each product³

1. Custom data reduction required. Use ADQ1600 Development Kit.
2. Custom data reduction required. Use SDR14 Development Kit.
3. Contact an SP Devices' sales representative for information about optional configurations.

Important Information

Teledyne Signal Processing Devices Sweden AB (Teledyne SP Devices) reserve the right to make corrections, modifications, enhancements, improvements, and other changes to its products and services at any time and to discontinue any product or service without notice. Customers should obtain the latest relevant information before placing orders and should verify that such information is current and complete. All products are sold subject to Teledyne SP Devices' general terms and conditions supplied at the time of order acknowledgment.

Teledyne SP Devices warrants that each product will be free of defects in materials and workmanship, and conform to specifications set forth in published data sheets, for a period of one (1) year. The warranty commences on the date the product is shipped by Teledyne SP Devices. Teledyne SP Devices' sole liability and responsibility under this warranty is to repair or replace any product which is returned to it by Buyer and which Teledyne SP Devices determines does not conform to the warranty. Product returned to Teledyne SP Devices for warranty service will be shipped to Teledyne SP Devices at Buyer's expense and will be returned to Buyer at Teledyne SP Devices' expense. Teledyne SP Devices will have no obligation under this warranty for any products which (i) has been improperly installed; (ii) has been used other than as recommended in Teledyne SP Devices' installation or operation instructions or specifications; or (iii) has been repaired, altered or modified by entities other than Teledyne SP Devices. The warranty of replacement products shall terminate with the warranty of the product. Buyer shall not return any products for any reason without the prior written authorization of Teledyne SP Devices.

In no event shall Teledyne SP Devices be liable for any damages arising out of or related to this document or the information contained in it.

TELEDYNE SP DEVICES' EXPRESS WARRANTY TO BUYER CONSTITUTES TELEDYNE SP DEVICES' SOLE LIABILITY AND THE BUYER'S SOLE REMEDY WITH RESPECT TO THE PRODUCTS AND IS IN LIEU OF ALL OTHER WARRANTIES, LIABILITIES AND REMEDIES. EXCEPT AS THUS PROVIDED, TELEDYNE SP DEVICES DISCLAIMS ALL WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR NON-INFRINGEMENT.

TELEDYNE SP DEVICES DOES NOT INDEMNIFY, NOR HOLD THE BUYER HARMLESS, AGAINST ANY LIABILITIES, LOSSES, DAMAGES AND EXPENSES (INCLUDING ATTORNEY'S FEES) RELATING TO ANY CLAIMS WHATSOEVER. IN NO EVENT SHALL TELEDYNE SP DEVICES BE LIABLE FOR SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING LOST PROFIT, LOST DATA AND THE LIKE, DUE TO ANY CAUSE WHATSOEVER. NO SUIT OR ACTION SHALL BE BROUGHT AGAINST TELEDYNE SP DEVICES MORE THAN ONE YEAR AFTER THE RELATED CAUSE OF ACTION HAS ACCRUED. IN NO EVENT SHALL THE ACCRUED TOTAL LIABILITY OF TELEDYNE SP DEVICES FROM ANY LAWSUIT, CLAIM, WARRANTY OR INDEMNITY EXCEED THE AGGREGATE SUM PAID TO SP BY BUYER UNDER THE ORDER THAT GIVES RISE TO SUCH LAWSUIT, CLAIM, WARRANTY OR INDEMNITY.

Worldwide Sales and Technical Support

www.spdevices.com

Teledyne SP Devices Corporate Headquarters

Teknikringen 6
SE-583 30 Linköping
Sweden

Phone: +46 (0)13 465 0600
Fax: +46 (0)13 991 3044
Email: info@spdevices.com

Copyright © 2017 Teledyne Signal Processing Devices Sweden AB. All rights reserved, including those to reproduce this publication or parts thereof in any form without permission in writing from Teledyne SP Devices.